

Ribbon Wetting Resin!

One of the first steps in using Ribbon is to wet the Ribbon with a plain unfilled enamel bonding resin. Most one bottle bonding resins contain organic solvents, dentin primers or self-etchants; these types of resins should not be used to wet the Ribbon. Use your favorite bonding system to bond to tooth structure but wet the Ribbon with an enamel bonding resin that does not contain solvents, primers and etchants. If your office does not have a suitable enamel bonding resin, try *Ribbon Wetting Resin*.

Ribbon® Wetting Resin

- * No alcohol, no organic solvents, primers, or acidic monomers
- * Unfilled resin for optimized fiber bonding
- * Works with all light and dual-cure resins
- * Can also be used as an instrument lubricant



Only \$28.00 for a 10 ml bottle

Ribbon products are sold directly by Ribbon, Inc.

800-624-4554 www.ribbon.com

Ref. 101410